

BECAUSE BREAKING UP IS HARD TO DO...



SPEX SamplePrep's Freezer/Mills® are cryogenic laboratory mills that chill samples in liquid nitrogen then pulverizes them with a magnetically driven impactor. They are recognized as the world's most effective laboratory mills.

TYPICAL SAMPLES

Plant and animal tissue, plastic and polymers, yeast, pharmaceuticals, food products, electronic components, textiles, hair & bones.

FEATURES AND BENEFITS

- • Sample vial maintains low level temperatures due to continuous immersion in liquid nitrogen.
- • Heat sensitive metabolites, isomers and complex molecules retain their composition as they are not degraded by heat or pressure.
- • Touch-screen interface for easy usage
- • Sample Vial reaches liquid nitrogen temperature in minutes
- • Maintenance-free design with only one moving part.
- • Backed by SPEX SamplePrep's 60 years of experience.



WATCH OUR VIDEOS!
Scan this QR code with your smartphone or visit www.youtube.com/sampleprep.

TOUCH SCREEN PROGRAM

Touch screen stores up to 20 user-defined grinding programs for easy quick recall. The new high resolution display is user-friendly and includes the following features:

- • Run history which is uploadable to a thumb drive via USB
- • Training videos and manuals that can be viewed on the screen.
- • A complete listing of accessories

6775 FREEZER/MILL®

The 6775 Freezer/Mill is a small cryogenic grinder with a self-contained liquid nitrogen tub and insulated case. Grinds samples from 0.1 to 5.0 grams using one standard vial or one microvial set.

6875 FREEZER/MILL®

The 6875 Freezer/Mill is a larger cryogenic grinder with a self-contained liquid nitrogen tub and insulated case. Optional auto-fill system is available to automatically fill the bath with liquid nitrogen. It grinds samples from 0.1 to 100 grams using one large vial, four small vials or four microvial sets.

FEATURES AVAILABLE IN BOTH UNITS

1. Safety features include LN sensor and lid interlock. If lid is opened during the grinding process grinding will stop.
2. Closed grinding vial with impactor prevents sample cross-contamination. Grinding vials made of alternative materials are also available.
3. Liquid nitrogen is required for operation. The sample's temperature is maintained because it NEVER leaves the enclosed liquid nitrogen bath.



6775
FREEZER/MILL®



6875
FREEZER/MILL®

To access application notes for the Freezer/Mill

Corporate Headquarters

SPEX SamplePrep
65 Liberty Street
Metuchen, NJ 08840
1-855-GET-SPEX

Tel: 732-623-0465
Fax: 732-906-2492
sampleprep@spex.com
www.spexsampleprep.com

Agence Nord:

Campus Effiscience
2 Rue Jean Perrin - Bât D
14460 Colombelles
tél : 02.31.34.50.74
fax : 02.31.34.55.17

Agence Sud:

Hôtel d'Entreprises de la Croix Rouge
Lot A4 - 10 Av de la Croix Rouge
84000 Avignon
tél : 04.90.27.17.95 fax : 04.90.27.17.52



Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél /fax : 03.88.04.01.81

www.deltalabo.fr
info@deltalabo.fr