



CG-500, CG-450 CG-400, CG-200 Freezer/Mill® Cryogenic Grinders



- Low maintenance design with no moving parts
- Equipped with liquid nitrogen level sensor, insulated case and lid safety interlock switch for operator protection (liquid nitrogen required for operation)
- Dual electromagnetic grinding technology
- Grinds samples in the 0.1 gram to 200 gram range
- Each grinding chamber holds one large grinding vial, one mid-size vial, four standard grinding vials, or four microvial sets

Some samples are impossible to grind at ambient temperatures. The Freezer/Mill is a unique grinding machine that uses liquid nitrogen to embrittle unconventional materials (that are difficult to grind at room temperature) and pulverizes them to a powder form without thermal degradation. Samples are placed in sealed cryogenic grinding vials then immersed in liquid nitrogen. The samples are cooled to cryogenic temperatures then ground by a magnetically shuttling impactor.

Since the vials are sealed, the cryogenic grinder maintains the integrity of its contents. Volatile compounds are retained, hazardous or critical samples are easily controlled, clean-up is simplified, and cross-sample contamination is eliminated.

Features



Intuitive color LCD touch screen controller



Magnetic grinding mechanism; no mechanical moving parts



Designed for vigorous horizontal milling



Continuous cooling with liquid nitrogen sample vial immersion

CG-500, CG-450, CG-400, CG-200 Freezer/Mill®

Typical Applications Include

DNA/RNA extraction, forensic/pharmaceutical testing, RoHS/WEEE, consumer product testing, food safety, materials research, and medical research.

Sample Types Include



Plastics & Polymers



Animal Tissue



Pharmaceuticals



Cannabis



Yeast


Before & After Samples – Bottle Caps



Specifications – CG-500 Freezer/Mill

High-capacity cryogenic grinder with dual grinding and cooling chambers. It can grind and cool up to 200 g of samples (0.5 to 100 g per chamber). Each grinding chamber holds one large grinding vial, one mid-size vial, four small grinding vials, or four microvial sets each. Comes equipped with liquid nitrogen auto-fill.

Voltage	115 VAC, 60 Hz or 230 VAC, 50 Hz
Dimensions	50.8 x 54.6 x 64.8 cm (20.0 x 21.5 x 25.5")
Net Weight	35.4 kg (78 lb)
Warranty	3 years




CG-500, CG-450, CG-400, CG-200 Freezer/Mill®

Specifications – CG-450 and CG-400 Freezer/Mill

Large cryogenic grinder that accommodates up to 100 g of samples. It has one dual-purpose chamber that pre-cools and grinds samples. Its grinding chamber holds one large grinding vial, one mid-size vial, four small grinding vials, or four microvial sets. The CG-450 comes equipped with liquid nitrogen auto-fill. The CG-400 is manually filled with liquid nitrogen.

Voltage	115 VAC, 60 Hz or 230 VAC, 50 Hz
Dimensions	52.0 x 55.0 x 46.0 cm (20.5 x 21.5 x 18.0")
Net Weight	20 kg (44 lb)
Warranty	3 years



Specifications – CG-200 Freezer/Mill

Compact cryogenic mill that accommodates samples from 0.1 to 5 g. It has one dual purpose chamber that pre-cools and grinds samples. It can hold one small grinding vial, or three microvials.

Voltage	115 VAC, 60 Hz or 230 VAC, 50 Hz
Dimensions	48.0 x 27.0 x 32.0 cm (19.0 x 10.5 x 12.8")
Net Weight	8.6 kg (18.9 lb)
Warranty	3 years



Need a Demo?

- Virtual sample processing – virtual demonstration of the equipment and sample processing
- Physical sample processing – physical testing of your sample with shared results
- Physical equipment testing – a demo unit can be sent to your facility to test before purchasing

Contact us today at sampleprep@antylia.com for details.

CG-500, CG-450, CG-400, CG-200 Freezer/Mill®

Ordering Information

Description	Item Number
CG-500 Freezer/Mill, Dual Chamber Cryogenic Grinder, 115 VAC, 60 Hz	04500-42
CG-500 Freezer/Mill, Dual Chamber Cryogenic Grinder, 230 VAC, 50 Hz	04577-94
CG-400 Freezer/Mill, High-Capacity Cryogenic Grinder, 115 VAC, 60 Hz	04500-40
CG-400 Freezer/Mill, High-Capacity Cryogenic Grinder, 230 VAC, 50 Hz	04577-92
CG-450 Freezer/Mill, High Capacity Cryogenic Grinder with Auto-Fill, 115 VAC, 60 Hz	04500-41
CG-450 Freezer/Mill, High Capacity Cryogenic Grinder with Auto-Fill, 230 VAC, 50 Hz	04577-93
CG-200 Freezer/Mill, Compact Cryogenic Grinder, 115 VAC, 60 Hz	04500-39
CG-200 Freezer/Mill, Compact Cryogenic Grinder, 230 VAC, 50 Hz	04577-77

Agence Nord:

ZA Object'Ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.detalabo.fr
info@detalabo.fr