

**DELTA LABO**

Plug N'Work – Campus EffiScience  
Bât D - 2 Rue Jean Perrin  
14460 COLOMBELLES  
Tél : 02.31.34.50.74 / Fax.: 02.31.34.55.17

**DELTA LABO**

Hôtel d' Entreprises de La Croix Rouge – Lot A4  
10 Avenue de La Croix Rouge  
84000 AVIGNON – FRANCE  
Tel. 04.90.27.17.95 / Fax.: 04.90.27.17.52

**10400-0 SETA SMOKE POINT APPARATUS**

**ASTM D1322; IP 57; BS 2000 Part 57;  
ISO-DIS 3014; NF M07-028**

Comprising a tripod stand, lamp body, and chimney. The measuring scale is viewed through a window on the front of the instrument.

**SUPPLIED WITH:** candle assembly and certified wick



Size (HxWxD):	47 x 21 x 19cm
Weight:	2.6kg

**ACCESSORIES**

REQUIRED	
10450-0	SETA MIRROR, for fixing to chimney to detect first indication of smoke.
OPTIONAL	
10410-0	SETA WICK TRIMMER, inserts wick and automatically sets to correct length.
10430-0	SETA INSTRUMENT CASE, for lamp assembly, stand and five candles.
CALIBRATION, VERIFICATION & TRAINING	
99850-0	MULTI-TEST VERIFICATION MATERIAL, ASTM D1322 Kerosine (for details see page 128).



**Support de mèche réf 10400-3**



**Verre de lampe avec bague de blocage réf 10400-003**



**Lot de 100 mèches, réf 10400-002**



**Bougie réf 10400-01**



**Caisse de transport réf 10430-2**