

## Flash Point Standards

Product Code	Test Method	Product Nature	Certified Value	Pack Size
CRM-COC-Low	COC ASTM D 92	Mineral Oil	85.0°C	3 x 80mL
CRM-COC-Mid	COC ASTM D 92	Mineral Oil	153.0°C	3 x 80mL
CRM-COC-High	COC ASTM D 92	Mineral Oil	262.0°C	3 x 80mL
CRM-PMCC-Low	PMCC ASTM D 93 Proc A	Mineral Oil	77.5°C	3 x 80mL
CRM-PMCC-Mid	PMCC ASTM D 93 Proc A	Mineral Oil	141.0°C	3 x 80mL
CRM-PMCC-High	PMCC ASTM D 93 Proc A	Mineral Oil	228.0°C	3 x 80mL
CRM-ABKR	Abel IP 170	Jet Aviation Fuel	41.9°C	250mL
CRM-FCLU	COC ASTM D 92	Lubricant	265°C	250mL
CRM-PMGO	PMCC ASTM D 93 Proc A	Diesel	58.2°C	250mL
CRM-PMLU	PMCC ASTM D 93 Proc B	Lubricant	203°C	250mL
CRM-TAKR	TAG ASTM D56	Jet Aviation Fuel	43.1°C	250mL

All values stated are nominal - 10/07/2011

	<b>Agence Nord :</b> Plug N'Work - Campus Effiscience 2 Rue Jean Perrin - Bât D 14460 Colombelles Tél: 02.31.34.50.74 Fax: 02.31.34.55.17
	<b>Agence Sud :</b> Hôtel d'Entreprises de La Croix Rouge - Lot A4 10 Av de la Croix Rouge - 84000 Avignon Tél : 04.90.27.17.95 Fax : 04.90.27.17.52
<a href="http://www.deltalabo.fr">www.deltalabo.fr</a>	