

**KDS100, KDS120**

<b>Syringe</b>	<b>Diameter</b>	<b>Minimum</b>	<b>Maximum</b>
10 µl	0.46 mm	0.1 µl/hr	126.5 µl/hr
25 µl	0.73 mm	0.1 µl/hr	318.8 µl/hr
50 µl	1.03 mm	0.2 µl/hr	625 µl/hr
100 µl	1.46 mm	1 µl/hr	1274 µl/hr
250 µl	2.3 mm	2 µl/hr	3164 µl/hr
500 µl	3.26 mm	3 µl/hr	6359 µl/hr
1 ml	4.61 mm	0.01 ml/hr	13.2 ml/hr
2.5 ml	7.28 mm	0.02 ml/hr	31.7 ml/hr
3 ml	8.59 mm	0.02 ml/hr	44.9 ml/hr
5 ml	10.3 mm	0.03 ml/hr	87.0 ml/hr
10 ml	14.57 mm	0.1 ml/hr	125.0 ml/hr
20 ml	19.05 mm	0.1 ml/hr	219.0 ml/hr
30 ml	21.59 mm	0.1 ml/hr	282.0 ml/hr
50 ml	28.9 mm	0.2 ml/hr	503.2 ml/hr
60 ml	26.6 mm	0.2 ml/hr	426.0 ml/hr

	<b>DELTA LABO</b> Créativa Bât A - Site Agroparc 200 Rue Michel de Montaigne BP 21221 84911 AVIGNON Cédex 9
	VOTRE PARTENAIRE LABORATOIRE

Tél / fax 04 90 22 04 27  
[www.deltalabo.fr](http://www.deltalabo.fr)