

Fuels Testing Standards

Product Code	Test Method	Product Nature	Certified Value	Pack Size
CRM-ADKR	Acidity ASTM D 3242	Kerosene	0.0021 mg KOH/g	250mL
CRM-APKR	Aniline Point ASTM D 611	Jet Aviation fuel	62.30°C	250mL
CRM-ACKR	Aromatics Content ASTM D 5580	Gasoline	28.25%	250mL
CRM-BEGA	Benzene Content ASTM D 4053 / D 5580	Gasoline	0.827 % Vol.	250mL
CRM-CNGO	Cetane Number ASTM D 613	Diesel	52.6	1000mL
CRM-FIKR	FIA Aromatics ASTM D 1319	Kerosene	15.18 % vol	250mL
CRM-OMGA	Motor Octane Number ASTM D 2700	Gasoline (Unleaded)	86.9	1000mL
CRM-ORGA	Res. Octane Number ASTM D 2699	Gasoline (Unleaded)	98.0	1000mL
CRM-VPGA	Reid Vapour Pres. ASTM D 4953A, 323A	Gasoline	82.6 kPa	250mL

All values stated are nominal - 12/07/2011

<p><u>Agence Nord:</u> ZA Object'ifs Sud - Lot A3 6 Allée Emilie du Châtelet 14123 lfs tél : 02.31.34.50.74 fax : 02.31.34.55.17</p>		<p><u>Agence Est:</u> Parc Club des Tanneries 2 Rue de la Faisanderie 67380 Lingolsheim tél : 03.88.04.01.81 fax : 03.68.93.01.52</p>
<p><u>Agence Sud:</u> Bât Le Venango. 392 Rue Jean Dausset AGROPARC - BP11575 84916 Avignon Cédex 9 tél : 04.90.27.17.95 fax : 04.90.27.17.52</p>	<p>www.deltalabo.fr info@deltalabo.fr</p>	