

Element / Chemical Standards

Product Code	Test Method	Product Nature	Certified Value	Pack Size
CRM-SUGOHI	Sulphur Content ASTM D 5433; ISO 2084	Diesel	309.0 mg/kg	50mL
CRM-SUGOLO	Sulphur Content ASTM D 5453; ISO 2084	Diesel	7.2 mg/kg	250mL
CRM-SUKR	ASTM D 3227-ISO 3012	Kerosene	5.18 mg/kg	50mL

All values stated are nominal - 12/07/2011

Agence Nord:

ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Sud:


Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

**www.deltalabo.fr
info@deltalabo.fr**

X-ray Standards



Agence Nord :
 Plug N'Work - Campus Effiscience
 2 Rue Jean Perrin - Bât D
 14460 Colombelles
 Tél: 02.31.34.50.74 Fax: 02.31.34.55.17

Agence Sud :
 Hôtel d'Entreprises de La Croix Rouge - Lot A4
 10 Av de la Croix Rouge - 84000 Avignon
 Tél : 04.90.27.17.95 Fax : 04.90.27.17.52

www.deltalabo.fr

Products				
Product Code	Test Method	Product Nature	Certified Value	Pack Size
PS-SRF-1		Gasoline	0-1000 mg/kg (0; 50; 150; 300; 500; 750;1000; mg/kg)	150mL
PS-SHPL-1		Diesel	0-700 mg/kg (0; 50; 100; 200; 350; 500; 700; mg/kg)	150mL
PS-SHP-2		Domestic Heating Oil	0-0,250 % m/m (0; 0.035; 0.050; 0.100; 0.150; 0.200; 0.250 % m/m)	150mL
PS-SHP-3		Petroleum	0-0,700 % m/m (0; 0.40; 0.070; 0.100; 0.400; 0.700 % m/m)	150mL
PS-SHP-4		Petroleum	0 - 2,50 % m/m (0; 0.400; 0.700; 1.00; 1.50; 2.00;2.50 %m/m)	150mL
PS-SHHP-5		Fuel Oil	0 - 3.00 % m/m (0; 0.500; 1.00; 1.50; 2.00;2.50 % m/m)	150mL
PS-PBRF-1		Gasoline (Leaded)	0 - 30mg/L (0; 5; 10; 15; 25; 30 mg/L)	150mL
PS-PBRF-2		Gasoline (Leaded)	0 - 0.20g/L (0; 0.1; 0.05; 0.10; 0.15; 0.20 g/L)	150mL