

## Distillation Standards

Product Code	Test Method	Product Nature	Certified Value	Pack Size
CRM-DIGO	Distillation ASTM D 86; NF EN ISO 3405	Diesel	164°C - 357°C	250mL
CRM-DIGA	Distillation ASTM D 86; NF EN ISO 3405	Gasoline (unleaded)	29.9°C - 184 °C	250mL
CRM-DIKR	Distillation ASTM D 86; NF EN ISO 3405	Jet Aviation Fuel	151°C - 269°C	250mL

All values stated are nominal - 12/07/2011

**Agence Nord:**

ZA Object'Ifs Sud - Lot A3  
6 Allée Emilie du Châtelet  
14123 Ifs  
tél : 02.31.34.50.74  
fax : 02.31.34.55.17

**Agence Sud:**

Bât Le Venango. 392 Rue Jean Dausset  
AGROPARC - BP11575  
84916 Avignon Cédex 9  
tél : 04.90.27.17.95 fax : 04.90.27.17.52



**Agence Est:**

Parc Club des Tanneries  
2 Rue de la Faisanderie  
67380 Lingolsheim  
tél : 03.88.04.01.81  
fax : 03.68.93.01.52

**[www.deltalabo.fr](http://www.deltalabo.fr)  
[info@deltalabo.fr](mailto:info@deltalabo.fr)**