

High Temperature Viscosity Standards



High Temperature Type S3S

Overview **Specification** **Shipping** **Related Products**

Paragon ISO 17025 / ISO Guide 34 certified high temperature viscosity standards are for the calibration and verification of glass capillary viscometers and other viscosity measuring equipment where operating temperature is controlled precisely.



No 0649



No 4589

Products	Kinematic mm ² /s (cSt)		Dynamic mPa.s (cP)			Density g/mL		SFS			
Product Code	20°C 68°F	25°C 77°F	37.78°C 100°F	40°C 104°F	50°C 122°F	60°C 140°F	80°C 176°F	98.89°C 210°F	100°C 212°F	135°C 275°F	150°C 302°F
S3S	4.596	4.059	3.058	2.926	2.423	2.047	1.530	1.212	1.197		0.7424
S6S	10.59	8.945	6.113	5.761	4.500	3.618	2.502	1.880	1.852		1.066
S20S	43.32	33.89	19.63	18.02	12.70	9.360	5.645	3.854	3.774		1.867
S60S	158.2	117.4	59.95	53.96	34.94	23.89	12.73	7.944	7.752		3.226
N100S	327.6	236.6	113.2	100.9	62.51	41.09	20.51	12.19	11.86		4.538
S200S	639.5	451.3	204.7	180.8	107.9	68.64	32.38	18.44	17.91		6.350
S600S	2141	1454	599.5	521.2	291.0	173.9	73.68	38.65	37.36		11.41
S2000S	7929	5214	1982	1699	889.6	499.7	188.8	90.27	86.83		22.16
S8000S	34390	22088	7883	6684	3319	1767	598.1	258.7	247.8		51.24
N30000S	128417	80886				5595					111.6
N600		1454				173.9				15.40	
N2000		5222				499.2				31.29	

Agence Nord:

ZA Object'Ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17

Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52



Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

www.deltalabo.fr
info@deltalabo.fr