

Dynamic mPa.s (cP)

Product	20.00°C 68.00°F	20.50°C 68.90°F	21.00°C 69.80°F	21.50°C 70.70°F	22.00°C 71.60°F	22.50°C 72.50°F	23.00°C 73.40°F	23.50°C 74.30°F	24.00°C 75.20°F	24.50°C 76.10°F	25.00°C 77.00°F
RTM1	0.3207	0.3191	0.3176	0.3160	0.3143	0.3128	0.3112	0.3097	0.3082	0.3066	0.3051
RTM2	0.5526	0.5493	0.5460	0.5427	0.5395	0.5362	0.5329	0.5297	0.5266	0.5235	0.5204
RTM3	1.028	1.020	1.012	1.005	0.9965	0.9890	0.9815	0.9740	0.9666	0.9599	0.9525
RTM4	2.345	2.319	2.294	2.269	2.245	2.221	2.197	2.174	2.151	2.129	2.107
RTM5	3.787	3.738	3.688	3.641	3.594	3.548	3.503	3.459	3.415	3.372	3.330
RTM6	5.507	5.424	5.341	5.261	5.182	5.105	5.030	4.956	4.884	4.813	4.745
RTM7	8.954	8.799	8.642	8.487	8.343	8.200	8.060	7.923	7.790	7.659	7.533
RTM8	12.55	12.31	12.06	11.83	11.61	11.39	11.18	10.97	10.76	10.56	10.37
RTM9	17.88	17.50	17.12	16.76	16.41	16.07	15.74	15.42	15.11	14.79	14.50
RTM10	25.37	24.81	24.26	23.73	23.21	22.71	22.22	21.75	21.29	20.83	20.40
RTM11	36.92	35.99	35.08	34.21	33.36	32.55	31.76	30.99	30.25	29.52	28.83
RTM12	73.11	71.06	69.08	67.17	65.33	63.55	61.83	60.16	58.56	57.02	55.52
RTM13	100.6	97.64	94.76	91.97	89.27	86.74	84.26	81.86	79.55	77.33	75.19
RTM14	136.5	132.3	128.3	124.4	120.6	117.1	113.6	110.3	107.0	103.9	100.9
RTM15	206.4	199.7	193.2	187.1	181.1	175.5	170.0	164.6	159.6	154.7	150.0
RTM16	283.6	274.2	265.1	256.3	247.9	239.9	232.1	224.7	217.5	210.7	204.1
RTM17	352.3	340.1	328.3	317.1	306.3	295.9	286.0	276.5	267.3	258.6	250.2
RTM18	423.6	408.9	394.8	381.3	368.3	355.9	343.9	332.4	321.5	310.9	300.8
RTM19	500.3	482.4	465.2	448.7	432.9	417.8	403.3	389.4	376.1	363.4	351.1
RTM20	556.0	536.1	517.1	499.0	481.5	464.8	448.8	433.4	418.7	404.6	391.0
RTM21	708.2	682.4	657.7	634.0	611.3	589.6	568.8	548.9	529.9	511.6	494.0
RTM22	1044	1005	967.7	931.6	897.3	864.2	832.7	802.6	773.8	746.2	719.7
RTM23	1105	1068	1033	999.3	966.8	936.1	905.4	877.2	849.1	822.2	796.5
RTM24	1477	1420	1365	1313	1264	1216	1170	1127	1085	1045	1008
RTM25	1859	1785	1716	1650	1586	1525	1467	1413	1359	1309	1260
RTM26	2260	2169	2083	1999	1921	1845	1774	1705	1639	1577	1517
RTM27	2998	2878	2763	2653	2549	2449	2353	2262	2175	2092	2013
RTM28	3847	3688	3538	3393	3256	3125	3001	2882	2769	2660	2556
RTM29	6116	5858	5612	5379	5156	4944	4742	4549	4366	4190	4024
RTM30	7000	6702	6418	6149	5892	5648	5415	5193	4982	4780	4588
RTM31	7645	7318	7008	6713	6431	6164	5909	5665	5435	5214	5003
RTM32	9227	8828	8455	8099	7758	7434	7125	6831	6552	6285	6030
RTM33	11784	11270	10790	10330	9889	9470	9069	8689	8330	7985	7657
RTM34	12467	11930	11410	10920	10450	10000	9583	9182	8794	8429	8081
RTM35	15852	15120	14430	13780	13160	12560	12000	11470	10970	10490	10030
RTM36	23690	22630	21620	20670	19750	18890	18070	17290	16550	15840	15166
RTM37	30593	29220	27920	26690	25520	24400	23330	22330	21370	20460	19588
RTM38	63422	60500	57730	55100	52600	50230	47980	45840	43810	41880	40049
RTM39	115164	109800	104700	99900	95310	90980	86840	82930	79210	75680	72328

Agence Nord:
ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Est:
Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

Agence Sud:
Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr