

APPARATUS FOR THE DETERMINATION OF LOSS ON HEATING OF ASPHALTIC COMPOUNDS ASTM D 6 - D 1754 - IP 145

The apparatus consists of a convection oven manufactured as described in the ASTM E 145 specification (type IB), suitable for temperatures up to 180°C and equipped with a rotating shelf inside. An electronic controlled motor with gear reducer mounted on the top of the oven rotates the shelf at 5.5 rpm.



- Enamel finished steel case, benchtop model.
- Stainless steel interior with plenum fan for air circulation. Internal dimensions 380 x 380 x 380 mm. High efficiency thermal insulation.
- Stainless steel heater.
- Hinged door fitted with double wall window 300 x 200 mm, for internal view.
- Two openings, one on the bottom and one on the top for ventilation.
- Anodized aluminium horizontal rotating shelf 300 mm diameter, capable of supporting nine ASTM D 6 or three ASTM D 1754 sample containers. The shelf is suspended by a vertical shaft, centered with the horizontal interior dimensions and vertically located in the centre of the oven.
- Electronic controlled motor, 1/4 HP with built-in gear reducer that provides to rotate the shelf at a speed of 5.5 rpm.
- Microprocessor controlled thermoregulator with PID action and built-in digital display, 0.1°C accuracy. Probe: Pt100 RTD. Working range: from ambient to 180°C ± 1°C.
- Safety device against overheating.
- Easy access control box containing all the electronics: anodized aluminium panel with English written indications.
- English written user manual with installation instructions.
- For 220 V/50 Hz connections: 2000 W power consumption.
- Dimensions: 900 x 700 x 900 mm approx. Weight: 110 Kg approx.
- CE marked.

AD0006-100 Apparatus

ACCESSORIES

AD0006-A00 ASTM D 6 container
AD1754-A00 ASTM D 1754 container
AD1754-A01 ASTM D 1754 rotating shelf
CAL001 PT100 simulator
CAL003 Official Certificate for Pt100 simulator

CONSUMABLES

TA013C-N00 ASTM 13C thermometer (+155°C/+170°C)

Specifications may vary without notice.

The apparatus includes the items listed aside the picture, accessories etc. should be purchased separately.

Agence Nord:

ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17

Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52



Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

www.deltalabo.fr
info@deltalabo.fr